

## Product datasheet for **RC200137L3V**

### **TMEM38B (NM\_018112) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	TMEM38B (NM_018112) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM38B
Synonyms:	bA219P18.1; C9orf87; D4Ertd89e; OI14; TRIC-B; TRICB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018112
ORF Size:	873 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200137).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018112.1</a>
RefSeq Size:	3558 bp
RefSeq ORF:	876 bp
Locus ID:	55151
UniProt ID:	<a href="#">Q9NVV0</a>
Cytogenetics:	9q31.2
Domains:	DUF714
Protein Families:	Transmembrane


[View online »](#)

**MW:** 32.5 kDa

**Gene Summary:** This gene encodes an intracellular monovalent cation channel that functions in maintenance of intracellular calcium release. Mutations in this gene may be associated with autosomal recessive osteogenesis. [provided by RefSeq, Oct 2012]