

Product datasheet for **RC200143A1V**

Human SSBP3 (NM_018070) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SSBP3 (NM_018070) AAV Particle
Tag:	Myc-DDK
Symbol:	SSBP3
Synonyms:	CSDP; SSDP; SSDP1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200143 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGGATCGCC**

ATGTTTGCCAAAGGCAAAGGCTCGGCCGTGCCCTCGGATGGGCAGGCTCGGGAAAAGTTAGCTTTATACG
TCTACGAATATTTACTGCACGTAGGAGCACAGAAATCTGCACAGACCTTCTTATCGGAGATTCGATGGGA
AAAAACATCACGTTGGGAGAACCCTGGGTTTTGCACTCGTGGTGGTGTGATTTTGGGACCTTTAC
TGTGCAGCTCCTGAAAGGAGAGACACTTGTGAACATTCAGTGAAGCAAAGCCTTTTCATGATTATAGT
CAGCAGCTGCCCGAGCCCGTGTGGCAACATTCACCAACGATGGGATGCCGGGAGGCCCATCCC
GCCAGGTTTTCTTCAGGGTCTCCGGGTCACAGCCCTCGCCGACGCACAGCCTCCACCTCACAATCCT
AGCAGCATGATGGGACCCACAGTCAGCCTCCGGGAGGAGTTCCTGGGACACAGCCATTGCTGCCAATT
CTATGGATCCACACGACAACAAGGCCACCCCAACATGGGAGGATCAATGCAGAGAATGAACCCTCCCCG
AGGCATGGGGCCCATGGGTCCCGGCCACAGAATTACGGCAGCGGCATGAGACCACCACCAACTCCCTC
GGCCCCGCATGCCCGGATTAACATGGGCCCGGAGCTGGCAGACCCTGGCCAATCCTAACAGTGCTA
ACTCAATTCATACTCCTCCTCATCACCTGGTACCTATGTGGGACCCCTGGTGGTGGCGGTCTCCAGG
AACACCATTATGCCAGTCCCGCAGATTCAACAAATTCAGTACAACATCTACACAATGATTAATCCA
GTGCCGCTGGAGGCAGCCGGTCCAACATCCCGATGGGTCCCGGCTCGGACGGTCCGATGGCGGCATGG
GTGGCATGGAGCCACCCACATGAATGGATCATTAGGGTCAGGCGACATAGACGGACTTCCAAAAATTC
TCTAACAAACATAAGTGGCATTAGCAATCCTCCAGGCACCCCTCGAGATGACGGCGAGCTAGGAGGGAAC
TTCTCCACTCCTTTCAGAACGACAATTATTCTCCAAGCATGACGATGAGTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence:	<p>>RC200143 protein sequence Red=Cloning site Green=Tags(s)</p> <p>MFAGKKGSAVPSDQGAREKLALVYVEYLLHVGAQKSAQTFLSEIRWEKNITLGEPPGFLHSWCVFWDLY CAAPERRDTCEHSSEAKAFHDYSAAAAPSPVLGNIPPNDGMPGGPIPPGFFQGPPGSQSPHAQPPPHNP SSMGPHSQPPGGVPGTQPLLNSMDPTRQQGHPNMGGSMQRMNPPRGMGMPGPGQNYGSGMRPPPNL GPAMPGINMGPAGRPWPNSANSIPYSSSSPGTYVGPFGGGGPPGTPIMPSPADSTNSSDNIYTMINP YPPGSSRSNFPMGPSDGPMMGGMEPHHMNGSLGSGDIDGLPKNSPNNISGISNPPGTPRDDGELGGN FLHSFQNDNYSPTSMTMSV</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Species:	Human
Serotype:	AAV-2
ACCN:	NM_018070
ORF Size:	1104 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_018070.3
RefSeq Size:	3199 bp
RefSeq ORF:	1107 bp
Locus ID:	23648
UniProt ID:	Q9BWW4
Cytogenetics:	1p32.3
MW:	38.2 kDa