

## Product datasheet for **RC200151L2V**

### UNC45A (NM\_018671) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	UNC45A (NM_018671) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UNC45A
Synonyms:	GC-UNC45; GCUNC-45; GCUNC45; IRO039700; OOHE; SMAP-1; SMAP1; UNC-45A
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_018671
ORF Size:	2832 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200151).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018671.3</a>
RefSeq Size:	3585 bp
RefSeq ORF:	2835 bp
Locus ID:	55898
UniProt ID:	<a href="#">Q9H3U1</a>
Cytogenetics:	15q26.1
MW:	102.9 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a regulatory component of the progesterone receptor/heat shock protein 90 chaperoning complex, which functions in the assembly and folding of the progesterone receptor. The encoded protein is thought to be essential for normal cell proliferation, and for the accumulation of myosin during development of muscle cells. [provided by RefSeq, Sep 2018]