

Product datasheet for **RC200155L3V**

PLEKHB2 (NM_001031706) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PLEKHB2 (NM_001031706) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PLEKHB2
Synonyms:	EVT2; EVT2, FLJ20783; FLJ20783; FLJ23679; FLJ30436; pleckstrin homology domain containing, family B (evectins) member 2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001031706
ORF Size:	927 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200155).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001031706.1 , NP_001026876.1
RefSeq Size:	1197 bp
RefSeq ORF:	929 bp
Locus ID:	55041
Cytogenetics:	2q21.1
MW:	34.2 kDa
Gene Summary:	Involved in retrograde transport of recycling endosomes.[UniProtKB/Swiss-Prot Function]



[View online »](#)