

Product datasheet for **RC200160A1V**

Human C9orf95 (NMRK1) (NM_017881) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C9orf95 (NMRK1) (NM_017881) AAV Particle
Tag: Myc-DDK
Symbol: C9orf95
Synonyms: bA235O14.2; C9orf95; NRK1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200160 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAAACATTTATCATTGGAATCAGTGGTGTGACAAACAGTGGCAAACAACACTGGCTAAGAATTTGC
 AGAAACACCTCCCAAATTGCAGTGTCTATCTCAGGATGATTTCTTCAAGCCAGAGTCTGAGATAGAGAC
 AGATAAAATGGATTTTTCAGTACGATGTGCTTGAAGCACTTAACATGGAAAAATGATGTCAGCCATT
 TCCTGCTGGATGGAAGCGCAAGACACTCTGTGGTATCAACAGACCAGGAAAGTCTGAGGAAATCCCA
 TTTAATCATCGAAGGTTTTCTCTTTTAATTATAAGCCCTTGACACTATATGAATAGAAGCTATT
 CCTGACGATTCCATATGAAGAATGTAAAGGAGGAGGTACAAGGGTCTATCAGCCTCCAGACTCTCCG
 GGATACTTTGATGGCCATGTGTGGCCCATGTATCTAAAGTACAGACAAGAAATGCAGGACATCACATGGG
 AAGTTGTGTACCTGGATGGAACAAAATCTGAAGAGGACCTCTTTTGAAGTATATGAAGATCTAATACA
 AGAACTAGCAAAGCAAAGTGTGCAAGTGACAGCA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200160 protein sequence
 Red=Cloning site Green=Tags(s)

MKTFIIGISGVTNSGKTTLAKNLQKHLNCSVISQDDFFKPESEIETDKNGFLQYDVLEALNMEKMMSAI
 SCWMESARHSVSTDQESAEEIPILIEGFLFNKPLDTIWNRSYFLTIPYEECKRRRSTRVYQPPDSP
 GYFDGHVWPMYLYRQEMQDITWEVVYLDGKSEEDLFLQVYEDLIQELAKQKCLQVTA

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_017881
ORF Size:	597 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_017881.1</u>
RefSeq Size:	1207 bp
RefSeq ORF:	600 bp
Locus ID:	54981
UniProt ID:	<u>Q9NWW6</u>
Cytogenetics:	9q21.13
MW:	23.2 kDa