

Product datasheet for **RC200160L4V**

C9orf95 (NMRK1) (NM_017881) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C9orf95 (NMRK1) (NM_017881) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C9orf95
Synonyms:	bA235O14.2; C9orf95; NRK1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017881
ORF Size:	597 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200160).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017881.1
RefSeq Size:	1207 bp
RefSeq ORF:	600 bp
Locus ID:	54981
UniProt ID:	Q9NWW6
Cytogenetics:	9q21.13
Protein Pathways:	Nicotinate and nicotinamide metabolism
MW:	23.2 kDa



[View online »](#)

Gene Summary:

Nicotinamide adenine dinucleotide (NAD⁺) is essential for life in all organisms, both as a coenzyme for oxidoreductases and as a source of ADP-ribosyl groups used in various reactions. Nicotinic acid and nicotinamide, collectively known as niacin, are the vitamin precursors of NAD⁺. Nicotinamide riboside kinases, such as NRK1, function to synthesize NAD⁺ through nicotinamide mononucleotide using nicotinamide riboside as the precursor (Bieganski and Brenner, 2004 [PubMed 15137942]).[supplied by OMIM, Mar 2008]