

Product datasheet for RC200162

WDR5 (NM_052821) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDR5 (NM_052821) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	WDR5
Synonyms:	BIG-3; CFAP89; SWD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC200162 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGACGGAGGAGAAGAAGCCCGAGACCGAGGCCGCCAGAGCACAGCCAACCCCTTCGTCATCCGCCA
CTCAGAGCAAGCCTACACCTGTGAAGCCAACTATGCTCTAAAGTTCACCCTTGTGGCCACACCAAAGC
AGTGTCTCCGTGAAATTCAGCCGAATGGAGAGTGGCTGGCAAGTTCATCTGCTGATAAACTTATTTAA
ATTTGGGGCGGTATGATGGGAAATTTGAGAAAACCATATCTGGTCACAAGCTGGGAATATCCGATGTAG
CCTGGTCGTGAGATTCTAACCTTCTGTTTCTGCCTCAGATGACAAAACCTTGAAGATATGGGACGTGAG
CTCGGGCAAGTGTCTGAAAACCCCTGAAGGGACACAGTAATTATGTCTTTTGTGCAACTTCAATCCCCAG
TCCAACCTTATTGCTCAGGATCCTTTGACGAAAGCGTGAGGATATGGGATGTGAAAACAGGGAAGTGCC
TCAAGACTTTGCCAGCTCACTCGGATCCAGTCTCGGCCGTTTCAATTTAATCGTGATGGATCCTTGATAGT
TTCAAGTAGTATGATGGTCTCTGTGCGATCTGGGACACCGCCTCGGGCCAGTGCCTGAAGACGCTCATC
GATGACGACAACCCCGGTGTCTTTTGTGAAGTTCTCCCGAACGGCAAATACATCCTGGCCGCCACGC
TGGACAACACTCTGAAGCTCTGGGACTACAGCAAGGGGAAGTGCCTGAAGACGTACACTGGCCACAAGAA
TGAGAAATACTGCATATTTGCCAATTTCTGTTACTGGTGGGAAGTGGATTGTGTCTGGCTCAGAGGAT
AACCTTGTTCATCTGGAACCTTCAGACGAAAGAGATTGTACAGAACTACAAGGCCACACAGATGTCG
TGATCTCAACAGCTTGTACCCCAACAGAAAACATCATCGCCTCTGCTGCGCTAGAAAATGACAAAACAAT
TAAACTGTGGAAGAGTGACTGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC200162 protein sequence
Red=Cloning site Green=Tags(s)

MATEEKKPETEAARAQPTPSSSATQSKPTPVKPNYALKFTLAGHTKAVSSVKFSPNGEWLASSADKLIK
 IWGAYDGKFEKTSI SGHKLGI SDVAWSSDNLV SASDDKTLKIWDVSSGKCLKTLKGHSNYVFCNFPQ
 SNLIVSGSFDES VRIWDVKTGKCLKTLP AHSDPVSAVHFNRDGLIVSSSYDGLCRIWDTASGQCLKTLI
 DDDNPPVSFVKFSPNGKYILAATLDNTLKLWDYSK GKCLKTYTGHKNEKYCIFANF SVTGGKWIVSGSED
 NLVYIWNLQTK EIVQKLQGH TDVVI STACHPTENIIASAALENDKTIKLWKS DC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6393_b11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_052821

ORF Size: 1002 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_052821.4
RefSeq Size:	3078 bp
RefSeq ORF:	1005 bp
Locus ID:	11091
UniProt ID:	P61964
Cytogenetics:	9q34.2
Protein Families:	Druggable Genome
MW:	36.6 kDa
Gene Summary:	<p>This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This protein contains 7 WD repeats. Alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008]</p>


Product images:



Circular map for RC200162



Western blot validation of overexpression lysate (Cat# [LY409441]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC200162 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified WDR5 protein (Cat# [TP300162]). The protein was produced from HEK293T cells transfected with WDR5 cDNA clone (Cat# RC200162) using MegaTran 2.0 (Cat# [TT210002]).