

## Product datasheet for **RC200174L4V**

### CHCHD3 (NM\_017812) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CHCHD3 (NM_017812) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CHCHD3
Synonyms:	Mic19; MICOS19; MINOS3; PPP1R22
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017812
ORF Size:	681 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200174).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017812.1</a>
RefSeq Size:	1622 bp
RefSeq ORF:	684 bp
Locus ID:	54927
UniProt ID:	<a href="#">Q9NX63</a>
Cytogenetics:	7q32.3-q33
MW:	26.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is an inner mitochondrial membrane scaffold protein. Absence of the encoded protein affects the structural integrity of mitochondrial cristae and leads to reductions in ATP production, cell growth, and oxygen consumption. This protein is part of the mitochondrial contact site and cristae organizing system (MICOS). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2015]