

## Product datasheet for **RC200181L2V**

### **FAM46A (TENT5A) (NM\_017633) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FAM46A (TENT5A) (NM_017633) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TENT5A
Synonyms:	C6orf37; FAM46A; OI18; XTP11
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_017633
ORF Size:	1326 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200181).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017633.2</a>
RefSeq Size:	5617 bp
RefSeq ORF:	1329 bp
Locus ID:	55603
UniProt ID:	<a href="#">Q96IP4</a>
Cytogenetics:	6q14.1
MW:	49.7 kDa



[View online »](#)

**Gene Summary:** Probable nucleotidyltransferase that may act as a non-canonical poly(A) RNA polymerase.  
[UniProtKB/Swiss-Prot Function]