

Product datasheet for **RC200199A1V**

Human DKK3 (NM_013253) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	DKK3
Synonyms:	REIC; RIG
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC200199 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGCGGCTTGGGGCCACCCTGCTGTGCCTGCTGCTGGCGGCGGCGGTCCCCACGGCCCCGCGCCCCG
CTCCGACGGCGACCTCGGCTCCAGTCAAGCCCGGCCCGGCTCTCAGCTACCCGACGAGGAGGCCACCCT
CAATGAGATGTTCCGCGAGGTTGAGGAACTGATGGAGGACACGCAGCACAATTCGCGCAGCGCGGTGGAA
GAGATGGAGGCAGAAGAAGCTGCTGCTAAAGCATCATCAGAAGTGAACCTGGCAAACCTACCTCCAGCT
ATCACAATGAGACCAACACAGACACGAAGGTTGGAAATAATACCATCCATGTGCACCGAGAAATTCACAA
GATAACCAACAACCAGACTGGACAAATGGTCTTTTCAGAGACAGTTATCACATCTGTGGGAGACGAAGAA
GGCAGAAGGAGCCACGAGTGCATCATCGACGAGGACTGTGGGCCACGATGTACTGCCAGTTTGCCAGCT
TCCAGTACACCTGCCAGCCATGCCGGGGCCAGAGGATGCTCTGCACCCGGGACAGTGAAGTGTGTGGAGA
CCAGCTGTGTGTCTGGGGTCACTGCACCAAAATGGCCACCAGGGGCAGCAATGGGACCATCTGTGACAAC
CAGAGGACTGCCAGCCGGGGCTGTGCTGTGCTTCCAGAGAGGCTGCTGTTCCCTGTGTGCACACCCC
TGCCCGTGGAGGGCGAGCTTTGCCATGACCCCGCCAGCCGGCTTCTGGACCTCATCACCTGGGAGCTAGA
GCCTGATGGAGCCTTGGACCGATGCCCTTGTGCCAGTGGCCTCCTCTGCCAGCCCCACAGCCACAGCCTG
GTGTATGTGTGCAAGCCGACCTTCGTGGGGAGCCGTGACCAAGATGGGGAGATCCTGCTGCCAGAGAGG
TCCCCGATGAGTATGAAGTTGGCAGCTTCATGGAGGAGGTGCGCCAGGAGCTGGAGGACCTGGAGAGGAG
CCTGACTGAAGAGATGGCGCTGGGGAGCCTGCGGCTGCCGCCGCTGCACTGCTGGGAGGGGAAGAGATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>RC200199 protein sequence Red=Cloning site Green=Tags(s) MQR LGATLLCLLLAAAVPTAPAPAPTATSAPVKPGPALSYQEEATLNEMFREVEELMEDTQHKLRSAVE EMEAEEAAKASSEVNLANLPPSYHNETNTDTKVGNNITIHVHREIHKITNNQTGQMVFSETVITSVGDEE GRRSHECIIDEDCGPSMYCQFASFQYTCQPCRGQRM LCTR DSECCGQ LCVWGHCTKMATRGSNGTICDN QRDCQPG LCCAFQRGLLFPVCTPLPVEGELCHDPASRLDLITWELEPDGALDRCP CASGLLCQPHSHSL VYVCKPTFVGSRDQDGEILLPREVPDEYEVGSFMEEVRQELEDLERSLTEEMALGEPAAAAAALLGGEEI TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_013253
ORF Size:	1050 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_013253.4
RefSeq Size:	2755 bp
RefSeq ORF:	1053 bp
Locus ID:	27122
UniProt ID:	Q9UBP4
Cytogenetics:	11p15.3
MW:	38.3 kDa