

Product datasheet for **RC200200L3V**

ACOT7 (NM_007274) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ACOT7 (NM_007274) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ACOT7
Synonyms:	ACH1; ACT; BACH; CTE-II; hBACH; LACH; LACH1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007274
ORF Size:	1110 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200200).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007274.3
RefSeq Size:	1814 bp
RefSeq ORF:	1113 bp
Locus ID:	11332
UniProt ID:	O00154
Cytogenetics:	1p36.31
Domains:	4HBT
Protein Pathways:	Biosynthesis of unsaturated fatty acids



[View online »](#)

MW: 40.3 kDa

Gene Summary: This gene encodes a member of the acyl coenzyme family. The encoded protein hydrolyzes the CoA thioester of palmitoyl-CoA and other long-chain fatty acids. Decreased expression of this gene may be associated with mesial temporal lobe epilepsy. Alternatively spliced transcript variants encoding distinct isoforms with different subcellular locations have been characterized. [provided by RefSeq, Jul 2008]