

Product datasheet for **RC200218L4V**

ST6GALNAC6 (NM_013443) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ST6GALNAC6 (NM_013443) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ST6GALNAC6
Synonyms:	SIAT7-F; SIAT7F; ST6GALNACVI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013443
ORF Size:	999 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200218).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013443.3
RefSeq Size:	2474 bp
RefSeq ORF:	1002 bp
Locus ID:	30815
UniProt ID:	Q969X2
Cytogenetics:	9q34.11
Domains:	Glyco_transf_29
Protein Families:	Transmembrane



[View online »](#)

Protein Pathways: Glycosphingolipid biosynthesis - ganglio series, Metabolic pathways

MW: 38.1 kDa

Gene Summary: ST6GALNAC6 belongs to a family of sialyltransferases that modify proteins and ceramides on the cell surface to alter cell-cell or cell-extracellular matrix interactions (Tsuchida et al., 2003 [PubMed 12668675]).[supplied by OMIM, Mar 2008]