

## Product datasheet for **RC200221L2V**

### CRABP2 (NM\_001878) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CRABP2 (NM_001878) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CRABP2
Synonyms:	CRABP-II; RBP6
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001878
ORF Size:	414 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200221).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001878.2</a>
RefSeq Size:	1088 bp
RefSeq ORF:	417 bp
Locus ID:	1382
UniProt ID:	<a href="#">P29373</a>
Cytogenetics:	1q23.1
Domains:	lipocalin
Protein Families:	Druggable Genome, Transcription Factors


[View online »](#)

**MW:** 15.7 kDa

**Gene Summary:** This gene encodes a member of the retinoic acid (RA, a form of vitamin A) binding protein family and lipocalin/cytosolic fatty-acid binding protein family. The protein is a cytosol-to-nuclear shuttling protein, which facilitates RA binding to its cognate receptor complex and transfer to the nucleus. It is involved in the retinoid signaling pathway, and is associated with increased circulating low-density lipoprotein cholesterol. Alternatively spliced transcript variants encoding the same protein have been found for this gene.[provided by RefSeq, Dec 2010]