

## Product datasheet for **RC200227L2V**

### TOLLIP (NM\_019009) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TOLLIP (NM_019009) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TOLLIP
Synonyms:	IL-1RAcPIP
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_019009
ORF Size:	822 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200227).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019009.2</a>
RefSeq Size:	3665 bp
RefSeq ORF:	825 bp
Locus ID:	54472
UniProt ID:	<a href="#">Q9H0E2</a>
Cytogenetics:	11p15.5
Domains:	C2, CUE
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Toll-like receptor signaling pathway

**MW:** 30.3 kDa

**Gene Summary:** This gene encodes a ubiquitin-binding protein that interacts with several Toll-like receptor (TLR) signaling cascade components. The encoded protein regulates inflammatory signaling and is involved in interleukin-1 receptor trafficking and in the turnover of IL1R-associated kinase. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016]