

Product datasheet for **RC200244A1V**

Human CLEC4E (NM_014358) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CLEC4E (NM_014358) AAV Particle
Tag: Myc-DDK
Symbol: CLEC4E
Synonyms: CLECSF9; MINCLE
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200244 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAATTCATCTAAATCATCTGAAACACAATGCACAGAGAGAGGATGCTTCTCTCCAAATGTTCTTAT
GGACTGTTGCTGGGATCCCCATCTATTTCTCAGTGCCTGTTTCATCACCAGATGTGTTGTGACATTTTCG
CATCTTTCAAACCTGTGATGAGAAAAAGTTTCAGCTACCTGAGAATTTACAGAGCTCTCCTGCTACAAT
TATGGATCAGGTTTCAGTCAAGAATTGTTGTCCATTGAAGTGGGAATATTTCAATCCAGCTGCTACTTCT
TTTCTACTGACACCATTTCTGGGCGTTAAGTTTAAAGAACTGCTCAGCCATGGGGGCTCACCTGGTGGT
TATCAACTCACAGGAGGAGCAGGAATTCCTTTCTACAAGAAACCTAAAATGAGAGAGTTTTTTATTGGA
CTGTCAGACCAGGTTGTCGAGGGTCAGTGGCAATGGGTGGACGGCACACCTTTGACAAAGTCTCTGAGCT
TCTGGGATGTAGGGGAGCCCAACAACATAGCTACCCTGGAGGACTGTGCCACCATGAGAGACTCTTCAA
CCCAAGGCAAAATTGGAATGATGTAACCTGTTTCTCAATTATTTTCGGATTTGTGAAATGGTAGGAATA
AATCCTTTGAACAAAGGAAAATCTCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200244 protein sequence
Red=Cloning site Green=Tags(s)

MNSSKSSETQCTERGCFFSSQMFLWTAVGIPILFLSACFITRCVVTFRIFQTCDEKKFQLPENFTELSCYN
YSGSVKNCCPLNWEYFQSSCYFFSTDITISWALSLKNCSAMGAHLVVINSQEEQEFLSYKKPKMREFFIG
LSDQVVEGQWQWVDGTPLTSLSFWDVGEPPNIATLEDCAITMRDSSNPRQNWNDVTCFLNYFRICEMVGI
NPLNKGKSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_014358
ORF Size:	657 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014358.2</u>
RefSeq Size:	2158 bp
RefSeq ORF:	660 bp
Locus ID:	26253
UniProt ID:	<u>Q9ULY5</u>
Cytogenetics:	12p13.31
MW:	25.1 kDa