

Product datasheet for **RC200250L1V**

EBNA1 binding protein 2 (EBNA1BP2) (NM_006824) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EBNA1 binding protein 2 (EBNA1BP2) (NM_006824) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EBNA1 binding protein 2
Synonyms:	EBP2; NOBP; P40
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_006824
ORF Size:	918 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200250).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006824.1
RefSeq Size:	1397 bp
RefSeq ORF:	921 bp
Locus ID:	10969
UniProt ID:	Q99848
Cytogenetics:	1p34.2
MW:	34.9 kDa



[View online »](#)

Gene Summary: Required for the processing of the 27S pre-rRNA.[UniProtKB/Swiss-Prot Function]