

Product datasheet for **RC200253L3V**

ZFPL1 (NM_006782) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZFPL1 (NM_006782) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZFPL1
Synonyms:	D11S750; MCG4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006782
ORF Size:	930 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200253).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006782.3
RefSeq Size:	1400 bp
RefSeq ORF:	933 bp
Locus ID:	7542
UniProt ID:	O95159
Cytogenetics:	11q13.1
Protein Families:	Transcription Factors, Transmembrane
MW:	34.1 kDa



[View online »](#)

Gene Summary:

Required for cis-Golgi integrity and efficient ER to Golgi transport. Involved in the maintenance of the integrity of the cis-Golgi, possibly via its interaction with GOLGA2/GM130. [UniProtKB/Swiss-Prot Function]