

Product datasheet for **RC200255A1V**

Human TMED1 (NM_006858) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TMED1 (NM_006858) AAV Particle
Tag: Myc-DDK
Symbol: TMED1
Synonyms: IL1rl1; IL1RL1LG; p24g1; Tp24
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200255 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATGGCGGCCGCGCGGCCCTAGCCCTGGCCTTGTGGCTACTAATGCCACCAAGTGGAGGTGGGAGGGG
CGGGGCCCCGCCAATCCAGGACGGTGAGTTCACGTTCTGTTGCCGCGGGGAGGAAGCAGTGTCTTA
CCAGTCCGCGCGGCCAACGCAAGCCTCGAGACCGAATACCAGGTGATCGGAGGTGCTGGACTGGACGTG
GACTTCACGCTGGAGAGCCCTCAGGGCGTGCTGTTGGTCAGCGAGTCCCGCAAGGCTGATGGGGTACACA
CGGTGGAGCCAAACGGAGGCCGGGACTACAAGCTGTGCTTTGACAACCTCTCAGCACCATCTCCGAGAA
GCTGGTGTCTTTGAACTGATCTTTGACAGCCTCCAGGATGACGAGGAGGTGCAAGGATGGGCAGAGGCT
GTGGAGCCCGAGGAGATGCTGGATGTTAAATGGAGGACATCAAGGAGTCCATTGAGACCATGCGGACCC
GGCTGGAGCGCAGCATCCAGATGCTCACGCTACTGCGGGCCTTCGAGGCACGTGACCGCAACCTGCAAGA
GGGCAACTTGGAGCGGTCAACTTCTGGTCAGCTGTCAACGTGGCGGTGCTGCTGCTGGTGGCTGTGCTG
CAGGTCTGCACGCTCAAGCGCTTCTCCAGGACAAGCGCCGGTGCCACG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200255 protein sequence
Red=Cloning site Green=Tags(s)

MMAAGALALALWLLMPPVEVGAGPPPIQDGEFTFLLPAGRKQCFYQSAPANASLETEYQVIGGAGLDV
DFTLESPQGVLLVSESRAKADGVHTVEPTEAGDYKLCFDFNSFTISEKL VFFELIFDSLQDDEEVEGWAEA
VEPEEMLDVKMEDIKESIETMRRLERSIQMLTLLRAFEARDRLQEGNLERVNFWSAVNAVLLLVAVL
QVCTLRFFQDKRPVPT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_006858
ORF Size:	681 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006858.2
RefSeq Size:	1739 bp
RefSeq ORF:	684 bp
Locus ID:	11018
UniProt ID:	Q13445
Cytogenetics:	19p13.2
MW:	25.2 kDa