

Product datasheet for **RC200262A1V**

Human SCAMP2 (NM_005697) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SCAMP2 (NM_005697) AAV Particle
Tag:	Myc-DDK
Symbol:	SCAMP2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200262 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCGGCTTTTCGACACCAACCCCTTCGCGGACCCAGTGGATGTAAACCCCTTCAGGATCCCTCTGTGA
CCCAGCTGACCAACGCCCCGAGGGCGGCTGGCGGAATTCACCCCTTCTCAGAGACAAATGCAGCGAC
AACAGTTCCTGTACCCAATCCCTGGGTCTCACAGCCAGCGGTTCTCCAGCCATCAGTGAACCAACC
CAGCCGACCCCCAGGCCGTGGTGTCTGCAGCCAGGCAGGCTGCTCCGGCAGCAGGAAGAACTGGACA
GGAAAGCTGCCGAGCTGGAACGCAAGGAGCGGGAGCTGCAGAACTGTAGCCAATTGCATGTGAGACA
GAACAATGGCCCCCTTGCCCTCGTGGTGCCCTGTGAAGCCCTGCTTCTATCAGGATTTCTCCACAGAG
ATCCCTGCCGACTACCGCGGATATGCAAGATGCTCTACTATCTGTGGATGTTGCATTCACTGACTCTGT
TTCTGAACCTGCTTGCCCTGCCTGGCCTGGTTCTCGGGCAACAGCTCCAAGGGAGTGGACTTTGGCCTCTC
CATCTGTGGTTTCTGATCTTCACTCCCTGTGCCTTCCTTTGTTGGTACCGACCCATCTATAAGGCCTTT
AGGTCCGACAACTCTTTCAGTCTCTTTGTGTTCTTCTTTGTATTTTTTGTCAAATAGGGATCTACATCA
TCCAGTTGGTTGGCATCCCTGGCCTGGGGGACAGCGGTTGGATTGCAGCCCTGTCTACACTGGATAATCA
TTCCCTGGCCATATCAGTCATCATGATGGTGGTGGCTGGCTTCTCACCTCTGTGCCGTGCTCTCAGTC
TTCCTCCTGCAGCGGGTGCACTCCCTCTACCGACGACAGGGGCCAGTTCAGCAGGCCCAGGAGGAGT
TTTCCCAGGGCATCTTCAGCAGCAGAACCTTCCACAGAGCTGCTTCATCTGCTGCCCAAGGAGCCTTCCA
GGGAAT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200262 protein sequence Red=Cloning site Green=Tags(s) MSAFDTNPFADPVDVNPFDPSVTQLTNAPQGGLAEFNPSETNAATTVPVTQLPGSSQPAVLQPSVEPT QPTPQAVVSAAQAGLLRQQEELDRKAAELERKERELQNTVANLHVRQNNWPPLPSWCPVKPCFYQDFSTE IPADYQRICKMLYYLWMLHSVTLFLNLLACLAWFSGNSSKGVDGFLSILWFLIFTPCAFLCWYRPIYKAF RSDNSFSFFVFFVFFCQIGIYIIQLVGIPGLGDSGWIAALSTLDNHSLAISVIMMVVAGFFTLCAVLSV FLLQRVHSLYRRTGASFQQAQEEFSQIGFSSRTFHRAASSAAQGAQGN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005697
ORF Size:	987 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_005697.3</u>
RefSeq Size:	1313 bp
RefSeq ORF:	990 bp
Locus ID:	10066
UniProt ID:	<u>O15127</u>
Cytogenetics:	15q24.1
MW:	36.6 kDa