

Product datasheet for **RC200268L3V**

Annexin A10 (ANXA10) (NM_007193) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Annexin A10 (ANXA10) (NM_007193) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ANXA10
Synonyms:	ANX14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007193
ORF Size:	972 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200268).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007193.3
RefSeq Size:	1447 bp
RefSeq ORF:	975 bp
Locus ID:	11199
UniProt ID:	Q9UJ72
Cytogenetics:	4q32.3
MW:	37.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined. [provided by RefSeq, Jul 2008]