

Product datasheet for **RC200276A1V**

Human TIP30 (HTATIP2) (NM_006410) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TIP30 (HTATIP2) (NM_006410) AAV Particle
Tag:	Myc-DDK
Symbol:	TIP30
Synonyms:	CC3; SDR44U1; TIP30
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200276 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCGGGCCTGCGGCGCTGAGCGCGGCGGCGGCGGCTGCTCTGGCGGCCGCCCTGCTCCTGCTGCGTC
GTGAGGACCCGGGGCGGGGGCTGGCCCCAGCATGGCCGAAACAGAAGCCCTGTCGAAGCTTCGGGAAGA
CTTCAGGATGCAGAATAAATCCGTCTTTATTTGGGCGCCAGCGGAGAAACCGGCAGAGTGCTCTTAAAG
GAAATCCTGGAGCAGGGCCTGTTTTCAAAGTCACGCTCATTGGCCGGAGGAAGCTCACCTTCGACGAGG
AAGCTTATAAAAAATGTGAATCAAGAAGTGGTGGACTTTGAAAAGTTGGATGACTACGCTCTGCCTTTCA
AGGTCATGATGTTGGATTCTGTTGCCTGGGTACCACCAGAGGGAAGCTGGGGCGGAGGGATTTGTTCTG
GTTGACCGAGATTATGTGCTGAAGTCTGCAGAGCTGGCAAAAGCTGGAGGGTGCAAACTTTCAACTTGC
TATCCTCTAAAGGAGCTGATAAATCAAGCAATTTTTATATCTACAAGTTAAGGGAGAAGTAGAAGCCAA
GGTTGAAGAATTAATAATTTGATCGTTACTCTGTATTTAGGCCTGGAGTTCTGTTATGTGATAGGCAAGAA
TCTCGCCAGGTGAATGGCTGGTTAGAAAGTCTTTGGCTCCTTACCAGACTCTTGGGCCAGGGGCATT
CTGTGCCTGTGGTGACCGTGGTTAGAGCAATGCTGAACAATGTGGTGAGACCAAGAGACAAGCAGATGGA
ACTGCTGGAGAACAAGGCCATCCATGACCTGGGAAAGCGCATGGCTCTCTCAAGCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200276 protein sequence Red=Cloning site Green=Tags(s) MAGPAALSAAAAAALAAALLLRREDPGPGAGPSMAETEALSKLREDFRMQNKSVFILGASGETGRVLLK EILEQGLFSKVTLIGRRKLTDEEAYKNVNQEVVDFEKLDDYASAFQGHVGFCCCLGTTTRGKAGAEFVR VDRDYVLKSAELAKAGGCKHFNLLSSKGADKSSNFLYLQVKGEVEAKVEELKFDRYSVFRPGVLLCDRQE SRPGEWLVRKFFGSLPDSWARGHSVPVTVVRAMLNNVVRPRDKQMELLENKAIHDLGKAHGSLKP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006410
ORF Size:	828 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006410.3, NP_006401.2</u>
RefSeq Size:	1404 bp
RefSeq ORF:	729 bp
Locus ID:	10553
UniProt ID:	<u>Q9BUP3</u>
Cytogenetics:	11p15.1
MW:	30.2 kDa