

Product datasheet for **RC200278L2V**

S100A1 (NM_006271) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	S100A1 (NM_006271) Human Tagged ORF Clone Lentiviral Particle
Symbol:	S100A1
Synonyms:	S100; S100-alpha; S100A
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006271
ORF Size:	282 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200278).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006271.1
RefSeq Size:	607 bp
RefSeq ORF:	285 bp
Locus ID:	6271
UniProt ID:	P23297
Cytogenetics:	1q21.3
MW:	10.5 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in stimulation of Ca²⁺-induced Ca²⁺ release, inhibition of microtubule assembly, and inhibition of protein kinase C-mediated phosphorylation. Reduced expression of this protein has been implicated in cardiomyopathies. [provided by RefSeq, Jul 2008]