

## Product datasheet for **RC200304L1V**

### **RAC2 (NM\_002872) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | RAC2 (NM_002872) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | RAC2   |
| Synonyms:                 | EN-7; Gx; HSPC022; IMD73A; IMD73B; IMD73C; p21-Rac2  |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-Myc-DDK (PS100064)  |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_002872  |
| ORF Size:                 | 576 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC200304).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_002872.4</a>  |
| RefSeq Size:              | 1538 bp  |
| RefSeq ORF:               | 579 bp   |
| Locus ID:                 | 5880   |
| UniProt ID:               | <a href="#">P15153</a>   |
| Cytogenetics:             | 22q13.1  |
| Domains:                  | ras, RAN, RAS, RHO, RAB  |
| Protein Families:         | Druggable Genome   |



[View online »](#)

**Protein Pathways:** Adherens junction, Axon guidance, B cell receptor signaling pathway, Chemokine signaling pathway, Colorectal cancer, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Leukocyte transendothelial migration, MAPK signaling pathway, Natural killer cell mediated cytotoxicity, Pancreatic cancer, Pathways in cancer, Regulation of actin cytoskeleton, VEGF signaling pathway, Viral myocarditis, Wnt signaling pathway

**MW:** 21.4 kDa

**Gene Summary:** This gene encodes a member of the Ras superfamily of small guanosine triphosphate (GTP)-metabolizing proteins. The encoded protein localizes to the plasma membrane, where it regulates diverse processes, such as secretion, phagocytosis, and cell polarization. Activity of this protein is also involved in the generation of reactive oxygen species. Mutations in this gene are associated with neutrophil immunodeficiency syndrome. There is a pseudogene for this gene on chromosome 6. [provided by RefSeq, Jul 2013]