

## Product datasheet for **RC200334L3V**

### FAM32A (NM\_014077) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM32A (NM_014077) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM32A
Synonyms:	OTAG-12; OTAG12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014077
ORF Size:	336 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200334).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014077.1</a>
RefSeq Size:	1464 bp
RefSeq ORF:	339 bp
Locus ID:	26017
UniProt ID:	<a href="#">Q9Y421</a>
Cytogenetics:	19p13.11
MW:	13.2 kDa


[View online »](#)

**Gene Summary:**

Isoform 1, but not isoform 2 or isoform 3, may induce G2 arrest and apoptosis. May also increase cell sensitivity to apoptotic stimuli.[UniProtKB/Swiss-Prot Function]