

## Product datasheet for **RC200343L4V**

### UXT (NM\_004182) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	UXT (NM_004182) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UXT
Synonyms:	ART-27; STAP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004182
ORF Size:	471 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200343).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004182.2</a>
RefSeq Size:	619 bp
RefSeq ORF:	474 bp
Locus ID:	8409
UniProt ID:	<a href="#">Q9UBK9</a>
Cytogenetics:	Xp11.23
Domains:	DUF232
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 18.2 kDa

**Gene Summary:** The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]