

## Product datasheet for **RC200351L3V**

### **RAB33A (NM\_004794) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RAB33A (NM_004794) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAB33A
Synonyms:	RabS10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004794
ORF Size:	711 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200351).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004794.2</a>
RefSeq Size:	1128 bp
RefSeq ORF:	714 bp
Locus ID:	9363
UniProt ID:	<a href="#">Q14088</a>
Cytogenetics:	Xq26.1
Protein Families:	Druggable Genome
MW:	26.6 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. It is GTP-binding protein and may be involved in vesicle transport. [provided by RefSeq, Jul 2008]