

Product datasheet for RC200352L4

RAB11A (NM_004663) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAB11A (NM_004663) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	RAB11A
Synonyms:	YL8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200352).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 8766
UniProt ID: [P62491](#)
Cytogenetics: 15q22.31
Domains: ras, RAN, RAS, RHO, RAB
Protein Families: Druggable Genome
Protein Pathways: Endocytosis
MW: 24.4 kDa

Gene Summary: The protein encoded by this gene belongs to the Rab family of the small GTPase superfamily. It is associated with both constitutive and regulated secretory pathways, and may be involved in protein transport. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2011]