

Product datasheet for **RC200361A1V**

Human TFIIF (GTF2F2) (NM_004128) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	TFIIF
Synonyms:	BTF4; RAP30; TF2F2; TFIIF
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC200361 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCGAGCGCGGGGAACCTCGACTTGACCGCGCCAAACAGAACACAGGAGTGTGGCTAGTCAAGGTTT
CTAAATATTTGTACAGCAATGGGCTAAAGCTTCTGGAAGAGGTGAAGTTGGGAACTCGCGATTGCCAA
GACTCAAGGAAGGACTGAGGTGTCATTTACTTTGAATGAGGATCTTGCAAATATTCATGATATTGGTGGA
AAACGAGCTTCAGTCAGTGCTCCTAGAGAATCCATTTGTCTTGCAAAGTGTGGAGGACAGACATTAA
CAGTATTTACTGAGAGCTCATCAGATAAGCTGTCATTGGAAGGAATAGTGGTACAAAGAGCTGAATGCCG
ACCAGCTGCCAGTGAAGTACATGCGATTAAAGATTGCAAATAGAAGAGTCTTCCAAACAGTGAGG
CTATCACAACAGCTGGACAAAGTTGTAACAACCAATTACAAACCTGTTGCTAATCATCAATACAATATCG
AATATGAAAGGAAAAAGAAAGACGGAAGCGAGCTCGAGCTGATAAACAACATGTTTTAGACATGCT
ATTTTCAGCCTTTGAGAAACATCAATACTATAATCTTAAGGACTTGGTGGACATCACAAAACACCTGTG
GTGTACCTGAAGGAAATCTTAAAGAAATTGGTGTTCAGAAATGTAAGGGATCCACAAAACACATGGG
AGCTGAAGCCAGGTACAGACACTATCAAGGAGAAGAAAGAGTGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200361 protein sequence
Red=Cloning site Green=Tags(s)

MAERGELDLTGAKQNTGVWLKVPKYLQQWAKASGRGEVGLRIAKTQGRTEVSFTLNEDLANIHDIGG
KPASVSAPREHPFVLQSVGGQTLTVFTESDDKLSLEGI VVQRAECRPAASENYMRKRLQIEESSKPV
LSQQLDKVVTNYKPVANHQYNI EYERKKKEDGKRARADKQHVLDM LFSAFEKHQYYNLKDLVDITKQPV
VYLKEILKEIGVQNVKGIHKN TWELKPEYRHYQGEKSD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Species:	Human
Serotype:	AAV-2
ACCN:	NM_004128
ORF Size:	747 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004128.1
RefSeq Size:	1469 bp
RefSeq ORF:	750 bp
Locus ID:	2963
UniProt ID:	P13984
Cytogenetics:	13q14.12-q14.13
MW:	28.4 kDa