

Product datasheet for **RC200368**

eIF2 alpha (EIF2S1) (NM_004094) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	eIF2 alpha (EIF2S1) (NM_004094) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	eIF2 alpha
Synonyms:	EIF-2; EIF-2A; EIF-2alpha; EIF2; EIF2A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC200368 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCGGTCTAAGTTGTAGATTTTATCAACACAAATTCCTGAGGTGGAAGATGTAGTGATGGTGAATG
TCAGATCCATTGCTGAAATGGGGCTTATGTCAGCTTGTGGAATACAACAACATTGAAGGCATGATTCT
TCTTAGTGAATTATCCAGAAGCGTATCCGTTCTATCAACAACTCATCCGAATTGGCAGGAATGAGTGT
GTGGTTGTCATTAGGGTGGACAAAGAAAAGGATATATTGATTTGTCAAAAAGAAGATTTCTCCAGAGG
AAGCAATCAAATGTGAAGACAAATTCACAAAATCCAAAAGTGTATAGCATTCTTCGTCATGTTGCTGA
GGTGTAGAAATACACCAAGGATGAGCAGCTGGAAGCCTATTCCAGAGGACTGCCTGGGTCTTTGATGAC
AAGTACAAGACCTGGATATGGTGCCTATGATGCATTTAAGCATGCAGTCTCAGACCCATCTATTTTGG
ATAGTTTGAATTTGAATGAAGATGAACGGGAAGTACTCATTAAATAATTAATAGGCGCTTGACCCACA
GGCTGTCAAAATTCGAGCAGATATTGAAGTGGCTTGTATGGTTATGAAGGCATTGATGCTGTAAGAA
GCCCTAAGAGCAGGTTTGAATGTTCTACAGAAAACATGCCATTAAGATTAATCTAATAGCTCCTCCTC
GGTATGTAATGACTACGACAACCCTGGAGAGAACAAGGCCCTTCTGTCTCAGTCAAGCTATGGCTGT
TATCAAAGAGAAGATTGAGGAAAAGAGGGGTGTGTTCAATGTTCAAATGGAGCCAAAGTGGTCACAGAT
ACAGATGAGACTGAACTTCCGAGGCAGATGGAGAGGCTTGAAAGAGAAAATGCCGAAGTGGATGGAGATG
ATGATGCAGAAGAAATGGAAGCCAAAGCTGAAGAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC200368 protein sequence
Red=Cloning site Green=Tags(s)

MPGLSCRFYQHKFPEVEDVVMNVRSIAEMGAYVSLLEYNNIEGMILLSELRRRIRINSINKLIRIGRNEC
 VVVIRVDKEKGYIDLKRRVSPPEAIKCEDKFTKSKTVYSILRHVAEVLEYTKDEQLES LFQRTAWVFDD
 KYKRPYGYAYDAFKHAVSDPSILDSLNLNEDEREVLINNINRRLTPQAVKIRADIEVACYGYEGIDAVKE
 ALRAGLNCSTENMPIKINLIAPPRYVMTTTTLERTEGLSVLSQAMAVIKEKIEEKRGVFNVMPEPKVVTD
 TDETELARQMERLERENAEVDGDDDAEEMEAKAED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6171_b01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_004094


ORF Size: 945 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_004094.5
RefSeq Size:	4165 bp
RefSeq ORF:	948 bp
Locus ID:	1965
UniProt ID:	P05198
Cytogenetics:	14q23.3
Domains:	S1
MW:	36.1 kDa
Gene Summary:	The translation initiation factor EIF2 catalyzes the first regulated step of protein synthesis initiation, promoting the binding of the initiator tRNA to 40S ribosomal subunits. Binding occurs as a ternary complex of methionyl-tRNA, EIF2, and GTP. EIF2 is composed of 3 nonidentical subunits, the 36-kD EIF2-alpha subunit (EIF2S1), the 38-kD EIF2-beta subunit (EIF2S2; MIM 603908), and the 52-kD EIF2-gamma subunit (EIF2S3; MIM 300161). The rate of formation of the ternary complex is modulated by the phosphorylation state of EIF2-alpha (Ernst et al., 1987 [PubMed 2948954]).[supplied by OMIM, Feb 2010]


Product images:




Circular map for RC200368



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY EIF2S1 (Cat# RC200368, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-EIF2S1 (Cat# [TA501374]). Positive lysates [LY401321] (100ug) and [LC401321] (20ug) can be purchased separately from OriGene.



Western blot validation of overexpression lysate (Cat# [LY401321]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC200368 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified EIF2S1 protein (Cat# [TP300368]). The protein was produced from HEK293T cells transfected with EIF2S1 cDNA clone (Cat# RC200368) using MegaTran 2.0 (Cat# [TT210002]).