

Product datasheet for **RC200368L4V**

eIF2 alpha (EIF2S1) (NM_004094) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	eIF2 alpha (EIF2S1) (NM_004094) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EIF2S1
Synonyms:	EIF-2; EIF-2A; EIF-2alpha; EIF2; EIF2A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004094
ORF Size:	945 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200368).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004094.4
RefSeq Size:	4165 bp
RefSeq ORF:	948 bp
Locus ID:	1965
UniProt ID:	P05198
Cytogenetics:	14q23.3
Domains:	S1
MW:	36.1 kDa



[View online »](#)

Gene Summary:

The translation initiation factor EIF2 catalyzes the first regulated step of protein synthesis initiation, promoting the binding of the initiator tRNA to 40S ribosomal subunits. Binding occurs as a ternary complex of methionyl-tRNA, EIF2, and GTP. EIF2 is composed of 3 nonidentical subunits, the 36-kD EIF2-alpha subunit (EIF2S1), the 38-kD EIF2-beta subunit (EIF2S2; MIM 603908), and the 52-kD EIF2-gamma subunit (EIF2S3; MIM 300161). The rate of formation of the ternary complex is modulated by the phosphorylation state of EIF2-alpha (Ernst et al., 1987 [PubMed 2948954]).[supplied by OMIM, Feb 2010]