

Product datasheet for **RC200370L2V**

FAM50A (NM_004699) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM50A (NM_004699) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM50A
Synonyms:	9F; DXS9928E; HXC-26; HXC26; MRXSA; XAP5
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004699
ORF Size:	1017 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200370).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004699.1
RefSeq Size:	1311 bp
RefSeq ORF:	1020 bp
Locus ID:	9130
UniProt ID:	Q14320
Cytogenetics:	Xq28
Domains:	XAP5
MW:	40.1 kDa



[View online »](#)

Gene Summary:

This gene belongs to the FAM50 family. The encoded protein is highly conserved in length and sequence across different species. It is a basic protein containing a nuclear localization signal, and may function as a DNA-binding protein or a transcriptional factor. [provided by RefSeq, Sep 2009]