

## Product datasheet for **RC200371L2V**

### DDX5 (NM\_004396) Human Tagged ORF Clone Lentiviral Particle

#### Product data:


Product Type:	Lentiviral Particles
Product Name:	DDX5 (NM_004396) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DDX5
Synonyms:	G17P1; HLR1; HUMP68; p68
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004396
ORF Size:	1842 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200371).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004396.2</a>
RefSeq Size:	3769 bp
RefSeq ORF:	1845 bp
Locus ID:	1655
UniProt ID:	<a href="#">P17844</a>
Cytogenetics:	17q23.3
Domains:	DEAD, helicase_C
Protein Pathways:	Spliceosome


[View online »](#)

**MW:** 69.1 kDa

**Gene Summary:** This gene encodes a member of the DEAD box family of RNA helicases that are involved in a variety of cellular processes as a result of its role as an adaptor molecule, promoting interactions with a large number of other factors. This protein is involved in pathways that include the alteration of RNA structures, plays a role as a coregulator of transcription, a regulator of splicing, and in the processing of small noncoding RNAs. Members of this family contain nine conserved motifs, including the conserved Asp-Glu-Ala-Asp (DEAD) motif, important to ATP binding and hydrolysis as well as RNA binding and unwinding activities. Dysregulation of this gene may play a role in cancer development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2017]

**Product images:**



[RC200371L2] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC200371L2V particle to overexpress human DDX5-mGFP fusion protein.