

Product datasheet for **RC200372L4V**

Drebrin (DBN1) (NM_004395) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Drebrin (DBN1) (NM_004395) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Drebrin
Synonyms:	D0S117E
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004395
ORF Size:	1947 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200372).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004395.2
RefSeq Size:	2942 bp
RefSeq ORF:	1950 bp
Locus ID:	1627
UniProt ID:	Q16643
Cytogenetics:	5q35.3
Domains:	ADF
MW:	71.4 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a cytoplasmic actin-binding protein thought to play a role in the process of neuronal growth. It is a member of the drebrin family of proteins that are developmentally regulated in the brain. A decrease in the amount of this protein in the brain has been implicated as a possible contributing factor in the pathogenesis of memory disturbance in Alzheimer's disease. At least two alternative splice variants encoding different protein isoforms have been described for this gene. [provided by RefSeq, Jul 2008]