

Product datasheet for **RC200406L4V**

Uridine Phosphorylase 1 (UPP1) (NM_003364) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Uridine Phosphorylase 1 (UPP1) (NM_003364) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UPP1
Synonyms:	UDRPASE; UP; UPASE; UPP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003364
ORF Size:	930 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200406).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003364.2
RefSeq Size:	1733 bp
RefSeq ORF:	933 bp
Locus ID:	7378
UniProt ID:	Q16831
Cytogenetics:	7p12.3
Domains:	PNP_UDP_1
Protein Pathways:	Drug metabolism - other enzymes, Metabolic pathways, Pyrimidine metabolism



[View online »](#)

MW: 33.9 kDa

Gene Summary: This gene encodes a uridine phosphorylase, an enzyme that catalyzes the reversible phosphorylation of uridine (or 2'- deoxyuridine) to uracil and ribose-1-phosphate (or deoxyribose-1-phosphate). The encoded enzyme functions in the degradation and salvage of pyrimidine ribonucleosides. [provided by RefSeq, Oct 2016]