

Species:	Human
Serotype:	AAV-2
ACCN:	NM_003973
ORF Size:	648 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003973.3 , NP_003964.2
RefSeq Size:	875 bp
RefSeq ORF:	648 bp
Locus ID:	9045
UniProt ID:	P50914
Cytogenetics:	3p22.1
MW:	23.5 kDa