

Product datasheet for **RC200463L1V**

Heme Oxygenase 1 (HMOX1) (NM_002133) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Heme Oxygenase 1 (HMOX1) (NM_002133) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Heme Oxygenase 1
Synonyms:	bK286B10; HMOX1D; HO-1; HSP32
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_002133
ORF Size:	864 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200463).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002133.1
RefSeq Size:	1606 bp
RefSeq ORF:	867 bp
Locus ID:	3162
UniProt ID:	P09601
Cytogenetics:	22q12.3
Domains:	Heme_oxygenase
Protein Families:	Druggable Genome, Transmembrane




[View online »](#)

Protein Pathways: Porphyrin and chlorophyll metabolism

MW: 32.8 kDa

Gene Summary: Heme oxygenase, an essential enzyme in heme catabolism, cleaves heme to form biliverdin, which is subsequently converted to bilirubin by biliverdin reductase, and carbon monoxide, a putative neurotransmitter. Heme oxygenase activity is induced by its substrate heme and by various nonheme substances. Heme oxygenase occurs as 2 isozymes, an inducible heme oxygenase-1 and a constitutive heme oxygenase-2. HMOX1 and HMOX2 belong to the heme oxygenase family. [provided by RefSeq, Jul 2008]

Product images:



[RC200463L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC200463L1V particle to overexpress human HMOX1-Myc-DDK fusion protein.