

Product datasheet for **RC200483L4V**

CRMP1 (NM_001313) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CRMP1 (NM_001313) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CRMP1
Synonyms:	CRMP-1; DPYSL1; DRP-1; DRP1; ULIP-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001313
ORF Size:	1716 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200483).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001313.3
RefSeq Size:	3074 bp
RefSeq ORF:	1719 bp
Locus ID:	1400
UniProt ID:	Q14194
Cytogenetics:	4p16.2
Domains:	Amidohydro_1
MW:	62.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]