

## Product datasheet for **RC200513L3V**

### **HAND1 (NM\_004821) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	HAND1 (NM_004821) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HAND1
Synonyms:	bHLHa27; eHand; Hxt; Thing1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004821
ORF Size:	645 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200513).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004821.1</a>
RefSeq Size:	1750 bp
RefSeq ORF:	648 bp
Locus ID:	9421
UniProt ID:	<a href="#">O96004</a>
Cytogenetics:	5q33.2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors
MW:	23.4 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to the basic helix-loop-helix family of transcription factors. This gene product is one of two closely related family members, the HAND proteins, which are asymmetrically expressed in the developing ventricular chambers and play an essential role in cardiac morphogenesis. Working in a complementary fashion, they function in the formation of the right ventricle and aortic arch arteries, implicating them as mediators of congenital heart disease. In addition, it has been suggested that this transcription factor may be required for early trophoblast differentiation. [provided by RefSeq, Jul 2008]