

## Product datasheet for **RC200526L4V**

### NDUFB10 (NM\_004548) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NDUFB10 (NM_004548) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NDUFB10
Synonyms:	PDSW
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004548
ORF Size:	516 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200526).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004548.1</a>
RefSeq Size:	741 bp
RefSeq ORF:	519 bp
Locus ID:	4716
UniProt ID:	<a href="#">O96000</a>
Cytogenetics:	16p13.3
Protein Pathways:	Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease



[View online »](#)

**MW:** 20.8 kDa

**Gene Summary:** Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.[UniProtKB/Swiss-Prot Function]