

Product datasheet for **RC200528A1V**

Human GAMT (NM_000156) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human GAMT (NM_000156) AAV Particle
Tag:	Myc-DDK
Symbol:	GAMT
Synonyms:	CCDS2; HEL-S-20; PIG2; TP53I2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200528 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGCGCCCCAGCGGACCCCCATCTTCGCGCCCGGCGAGAACTGCAGCCCCGCGTGGGGGGCGGCGC
CCGCGGCCTACGACGCAGCGGACACGCACCTGCGCATCCTGGGCAAGCCGGTGATGGAGCGCTGGGAGAC
CCCCTATATGCACGCGCTGGCCGCCGCCCTCCTCAAAGGGGGCGGGTCTGGAGGTGGGCTTTGGC
ATGGCCATCGCAGCGTCAAAGGTGCAGGAGCGCCATTGATGAGCATTGGATCATCGAGTGAATGACG
GCGTCTTCCAGCGGCTCCGGGACTGGGCCCCACGCGAGACACAAAGGTATCCCCTTGAAAGGCCTGTG
GGAGGATGTGGCACCCACCTGCCTGACGGTCACCTTTGATGGGATCCTGTACGACACGTACCCACTCTCG
GAGGAGACCTGGCACACACACCAGTTCAACTTCATCAAGAACCACGCCTTTCGCCTGCTGAAGCCGGGGG
GCGTCCTCACCTACTGCAACCTCACCTCCTGGGGGAGCTGATGAAGTCCAAGTACTCAGACATACCAT
CATGTTTGAGGAGACGCAGGTGCCCGCGCTGCTGGAGGCCGCTTCCGGAGGAGAACATCCGTACGGAG
GTGATGGCGCTGGTCCCACCGGCCGACTGCCGCTACTACGCCTTCCCACAGATGATCACGCCCTGGTGA
CCAAAGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200528 protein sequence Red=Cloning site Green=Tags(s) MSAPSATPIFAPGENCSPAWGAAPAAAYDAADTHLRILGKPVMERWETPYMHALAAAASSKGGRVLEVGFG MAIAASKVQEAPIDEHWIIECNDGVFQRLRDWAPRQTHKVIPLKGLWEDVAPTLPDGHFDGILYDTYPLS EETWHTHQFNFIKNHAFRLKPGGVLTYCNLTSGELMKSYSIDITIMFEETQVPALLEAGFRRENIRTE VMALVPPADCRYAFAFPQMITPLVTKG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_000156
ORF Size:	708 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_000156.5</u>
RefSeq Size:	1138 bp
RefSeq ORF:	711 bp
Locus ID:	2593
UniProt ID:	<u>Q14353</u>
Cytogenetics:	19p13.3
MW:	26.3 kDa