

Product datasheet for **RC200536L3V**

HEXIM1 (NM_006460) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HEXIM1 (NM_006460) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HEXIM1
Synonyms:	CLP1; EDG1; HIS1; MAQ1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006460
ORF Size:	1077 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200536).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006460.2
RefSeq Size:	4785 bp
RefSeq ORF:	1080 bp
Locus ID:	10614
UniProt ID:	O94992
Cytogenetics:	17q21.31
Protein Families:	Transcription Factors
MW:	40.6 kDa



[View online »](#)

Gene Summary: Expression of this gene is induced by hexamethylene-bis-acetamide in vascular smooth muscle cells. This gene has no introns. [provided by RefSeq, Jul 2008]