

Product datasheet for **RC200537A1V**

Human Stanniocalcin 2 (STC2) (NM_003714) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Stanniocalcin 2 (STC2) (NM_003714) AAV Particle
Tag:	Myc-DDK
Symbol:	Stanniocalcin 2
Synonyms:	STC-2; STCRP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200537 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGTGCCGAGCGGCTGGGCCAGTTCATGACCCTGGCTTTGGTGTGGCCACCTTTGACCCGGCGCGGG
GGACCGACGCCACCAACCCACCCGAGGGTCCCCAAGACAGGAGCTCCAGCAGAAAGGCCGCTGTCCCT
GCAGAATACAGCGGAGATCCAGCACTGTTTGGTCAACGCTGGCGATGTGGGTGTGGCGTGTGTAATGT
TTCGAGAACAACCTCTTGAGATTGGGGCTTACATGGGATTTGCATGACTTTTCTGCACAACGCTGGAA
AATTTGATGCCAGGGCAAGTCATTTCATCAAGACGCCCTTGAATGTAAGGCCACGCTCTGCGGCACAG
GTTTCGGCTGCATAAGCCGGAAGTGCCCGGCCATCAGGGAAATGGTGTCCAGTTGCAGCGGGAATGCTAC
CTCAAGCACGACCTGTGCGCGGCTGCCAGGAGAACACCCGGGTGATAGTGGAGATGATCCATTTCAAGG
ACTTGCTGCTGCACGAACCTACGTGGACCTCGTGAACCTGCTGCTGACCTGTGGGGAGGAGGTGAAGGA
GGCCATCACCCACAGCGTGCAGGTTCAAGTGTGAGCAGAACTGGGAAGCCTGTGCTCCATCTTGAGCTTC
TGACCTCGGCCATCCAGAAGCCTCCACGGCGCCCCCGAGCGCCAGCCCCAGGTGGACAGAACCAAGC
TCTCCAGGGCCCACCACGGGAAGCAGGACATCACCTCCAGAGCCAGCAGTAGGGAGACTGGCCGAGG
TGCCAAGGGTGAGCGAGGTAGCAAGAGCCACCAAACGCCCATGCCGAGGCAGAGTCGGGGCCTTGGG
GCTCAGGGACCTTCGGAAGCAGCGAGTGGGAAGACGAACAGTCTGAGTATTCTGATATCCGGAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200537 protein sequence Red=Cloning site Green=Tags(s) MCAERLGQFMTLALVLATFDPARGTDATNPPEGPQDRSSQKGRSLQNTAEIQHCLVNAGDVGCGVFEC FENNSCIEIRGLHGICMTFLHNAGKFDAQGKSF IKDALKCKAHALRHRFGCISRKCPAIREMVSQLQRECY LKHDLCAAAQENTRVIVEMIHFKDLLHEPYVDLVNLLTCGEEVKEAITHSVQVQCEQNWGSLCSILSF CTSAIQKPPTAPPERQPQVDRTKLSRAHHGEAGHHLPEPSSRETGRGAKGERGSKSHPNAHARGRVGGLG AQGPSGSSEWEDEQSEYSDIRR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003714
ORF Size:	906 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003714.2 , NP_003705.1
RefSeq Size:	5361 bp
RefSeq ORF:	909 bp
Locus ID:	8614
UniProt ID:	O76061
Cytogenetics:	5q35.2
MW:	33.2 kDa