

Product datasheet for **RC200549**

MSK1 (RPS6KA5) (NM_182398) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	MSK1
Synonyms:	MSK1; MSPK1; RLPK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC200549 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGGAGGAGGTGGCAGCAGCGCGCGCCGCGGGACCAGCGCGGACGGCGGCACGGAGGAGAGC
 AGCTCCTCACTGTCAAGCACGAGCTGCGGACTGCTAATTTGACAGGACATGCTGAGAAGGTGGGAATAGA
 AAATTTTGAGCTCCTGAAGGTCCTAGGAACTGGAGCTTATGAAAAAGTATTTCTAGTTCTGTAATAAAGT
 GGCCATGATACTGAAAGCTGTATGCCATGAAAGTTTAAAAAGGCAACAATCGTTCAAAGGCCAAAA
 CCACAGAGCATACAAGGACAGAACGACAAGTCTGGAACACATTAGGCAGTCGCCATTTTGGTAACATT
 ACATTATGCTTTCCAGACAGAAACCAACTTCATCTATTTTAGATTATATAAATGGTGGTGAACTTTTT
 ACTCATCTTTCTCAAAGAGAGCGTTTCACAGAGCATGAGGTGCAGATTTATGTTGGAGAGATTGTGCTTG
 CCCTCGAACATCTCCACAAGTTGGGATTATATATCGTGATTAAGCTTGAGAAATTCTACTTGATTC
 TAATGGCCATGTGGTCTGACAGATTTTGGTCTGAGTAAGGAGTTTGTGGCTGATGAAACTGAAAGAGCA
 TATTCCTTTGTGGAAGTATTGAATACATGGCACCAGATATTGTCAGAGGGGAGATTTCAGGACATGACA
 AGGCAGTTGACTGGTGGAGTTTGGGTGTTCTAATGTATGAATTACTAACTGGAGCATCTCCTTTCACTGT
 TGATGGAGAAAAAATCCCAAGCTGAGATATCTAGGAGAATATTAAGAGTGAGCCTCCATATCCCAA
 GAAATGAGTGCTTAGCGAAAGACCTAATTCAGCGTCTTTTGATGAAAGATCCCAAGAAGAGATTGGGAT
 GTGGTCCACGTGATGCAGATGAAATCAAAGAATCTCTTCTTTGAGAAAAATAAATGGGATGATTTAGC
 CGCCAAAAAGTGCTGCACCATTAAGCCAGTCATTCGAGATGAATTAGATGTGAGTAACCTTTCAGAA
 GAGTTCACAGAAATGGATCCCACTTATTCTCCCGCAGCCCTGCCCCAGAGTTCTGAGAAGCTGTTTCAGG
 GCTATTCCTTTGTGCTCCTTCCATCTATTCAAGCGTAATGCAGCTGTCATAGACCTCTTCAGTTTCA
 CATGGGAGTTGAACGTCCTGGAGTGACAAATGTTGCCAGGAGTGAATGATGAAGGACTCTCCATTCTAT
 CAACACTATGACCTAGATTTGAAGGACAAACCCCTGGGAGAAGGTAGTTTTTCAATTTGTCGAAAGTG
 TGCATAAAAAAAGTAACCAAGCTTTTGCAGTCAAATAATCAGCAAAAGGATGGAAGCCAATACTCAAAA
 GGAAATAACAGCTCTGAACTCTGTGAAGGACACCCCAATATTGTGAAGTTGCATGAAGTTTTTCATGAT
 CAGCTTCACACGTTTCTAGTGATGGAACCTCTGAATGGAGGAGAACTGTTTGAGCGCATTAAGAAAAAGA
 AGCACTTCAGTGAGACGGAAGCCAGCTACATCATGAGGAAGCTTGTTCAGCTGTAGCCACATGCATGA
 TGTTGGAGTGGTGCACAGGGATCTGAAACCTGAGGTA

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC200549 protein sequence
 Red=Cloning site Green=Tags(s)

MEEEGSSGGAAGTSADGGDGGEQLLTVKHELRTANLTGHAKEVGIENFELLKVLGTGAYGKVFLVRKIS
 GHDTGKLYAMKVLKKATIVQAKTTEHTRTERQVLEHIRQSPFLVTLHYAFQTETKLHLILDYINGGELF
 THLSQRERFTEHEVQIYVGEIVLALEHLHKLGIYRDIKENILLDSNGHVVLTDGFLSKEFVADETERA
 YSFCGTIEYMAPDIVRGDSGHDKAVDWWSLGVLMYELLTGASPFVTDGEKNSQAEISRRILKSEPPYPQ
 EMSALAKDLIQRLMKDPKKRLGCGPRDAEIKEHLFFQKINWDDLAAKKVPAPFKPVIDELDVSNFAE
 EFTEMDPTYSAPALPQSSEKLFQGSFVAPSILFKRNAVIDPLQFHMVERPGVTNVARSAMMKDSPFY
 QHYDLDLKDKPLGEGSFICRKCVMHKSNAFAVKIISKMEANTQKEITALKLCCEHPNIVKLHEVFHD
 QLHTFLVMELNNGGELFERIKKKHFFSETEASYIMRKLVSASVSHMHDVGVVHRDLKPEV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6153_b06.zip

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_182398

ORF Size: 1647 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_182398.2](#)


RefSeq Size: 2343 bp

RefSeq ORF: 1650 bp


Locus ID: 9252

UniProt ID:	<u>O75582</u>
Cytogenetics:	14q32.11
Protein Families:	Druggable Genome, Protein Kinase, Transcription Factors
Protein Pathways:	Bladder cancer, MAPK signaling pathway, Neurotrophin signaling pathway
MW:	61.8 kDa
Gene Summary:	Serine/threonine-protein kinase that is required for the mitogen or stress-induced phosphorylation of the transcription factors CREB1 and ATF1 and for the regulation of the transcription factors RELA, STAT3 and ETV1/ER81, and that contributes to gene activation by histone phosphorylation and functions in the regulation of inflammatory genes (PubMed:11909979, PubMed:12569367, PubMed:12763138, PubMed:9687510, PubMed:18511904, PubMed:9873047). Phosphorylates CREB1 and ATF1 in response to mitogenic or stress stimuli such as UV-C irradiation, epidermal growth factor (EGF) and anisomycin (PubMed:11909979, PubMed:9873047). Plays an essential role in the control of RELA transcriptional activity in response to TNF and upon glucocorticoid, associates in the cytoplasm with the glucocorticoid receptor NR3C1 and contributes to RELA inhibition and repression of inflammatory gene expression (PubMed:12628924, PubMed:18511904). In skeletal myoblasts is required for phosphorylation of RELA at 'Ser-276' during oxidative stress (PubMed:12628924). In erythropoietin-stimulated cells, is necessary for the 'Ser-727' phosphorylation of STAT3 and regulation of its transcriptional potential (PubMed:12763138). Phosphorylates ETV1/ER81 at 'Ser-191' and 'Ser-216', and thereby regulates its ability to stimulate transcription, which may be important during development and breast tumor formation (PubMed:12569367). Directly represses transcription via phosphorylation of 'Ser-1' of histone H2A (PubMed:15010469). Phosphorylates 'Ser-10' of histone H3 in response to mitogenics, stress stimuli and EGF, which results in the transcriptional activation of several immediate early genes, including proto-oncogenes c-fos/FOS and c-jun/JUN (PubMed:12773393). May also phosphorylate 'Ser-28' of histone H3 (PubMed:12773393). Mediates the mitogen- and stress-induced phosphorylation of high mobility group protein 1 (HMG1/HMG14) (PubMed:12773393). In lipopolysaccharide-stimulated primary macrophages, acts downstream of the Toll-like receptor TLR4 to limit the production of pro-inflammatory cytokines (By similarity). Functions probably by inducing transcription of the MAP kinase phosphatase DUSP1 and the anti-inflammatory cytokine interleukin 10 (IL10), via CREB1 and ATF1 transcription factors (By similarity). Plays a role in neuronal cell death by mediating the downstream effects of excitotoxic injury (By similarity). Phosphorylates TRIM7 at 'Ser-107' in response to growth factor signaling via the MEK/ERK pathway, thereby stimulating its ubiquitin ligase activity (PubMed:25851810).[UniProtKB/Swiss-Prot Function]


Product images:



Circular map for RC200549



Western blot validation of overexpression lysate (Cat# [LY403633]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC200549 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified RPS6KA5 protein (Cat# [TP300549]). The protein was produced from HEK293T cells transfected with RPS6KA5 cDNA clone (Cat# RC200549) using MegaTran 2.0 (Cat# [TT210002]).