

## Product datasheet for **RC200558L4V**

### TUSC2 (NM\_007275) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TUSC2 (NM_007275) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TUSC2
Synonyms:	C3orf11; FUS1; PAP; PDAP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007275
ORF Size:	330 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200558).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007275.1</a>
RefSeq Size:	1705 bp
RefSeq ORF:	333 bp
Locus ID:	11334
UniProt ID:	<a href="#">O75896</a>
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome
MW:	12.1 kDa


[View online »](#)

**Gene Summary:**

This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available. [provided by RefSeq, Jul 2008]