

Product datasheet for **RC200571A1V**

Human MPDU1 (NM_004870) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human MPDU1 (NM_004870) AAV Particle
Tag:	Myc-DDK
Symbol:	MPDU1
Synonyms:	CDGIF; HBEBP2BPA; Lec35; My008; PP3958; PQLC5; SL15; SLC66A5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200571 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCCGAGGCGGACGGACCGCTTAAACGGCTGCTCGTGCCGATTCTTTACCTGAGAAATGCTACG
ACCAACTTTTCGTTCACTGGGACTTGCTTCACGTCCCTGCCTCAAGATTCTCCTCAGCAAAGGCCTGGG
GCTGGGCATTGTGGCTGGCTCACTTCTAGTAAAGCTGCCCCAGGTGTTTAAATCCTGGGAGCCAAGAGT
GCTGAAGGTTGAGTCTCCAGTCTGTAATGCTGGAGCTAGTGGCATTGACTGGGACCATGGTCTACAGCA
TCACTAACAACTTCCATTCACTCTTGGGGTGAAGCCTTATTCCTGATGCTCCAGACGATCACCATCTG
CTTCTGCTCATGCACTACAGAGGACAGACTGTGAAAGGTGTCGCTTTCTCGCTTGCTACGGCCTGGTC
CTGCTGGTGTCTCTCACCTCTGACGCCCTTGACTGTAGTCAACCCTGCTCCAGGCCTCAATGTGCCTG
CTGTGGTGGTGGGAGGCTTCTCCAGGCAGCCACCACTACCACAACGGGCACACAGGCCAGCTCTCAGC
CATCACAGTCTTCTGCTGTTTGGGGCTCCCTGGCCGAATCTTCACTTCCATTAGGAAACCGGAGAT
CCCCTGATGGCTGGGACCTTTGTGCTCTCTCTGCAACGGCCTCATCGCCACCCAGCTGCTCTTCT
ACTGGAATGCAAAGCCTCCCCACAAGCAGAAAAAGGCGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200571 protein sequence Red=Cloning site Green=Tags(s) MAAEADGPLKRLLVPILLPEKCYDQLFVQWDLHVPCLKILLSKGLGLGIVAGSLLVKLPQVFKILGAKS AEGLSLQSVMLELVALTGTMVYSITNFPFSSWGEALFLMLQTITICFLVMHYRGQTVKGVAFLACYGLV LLVLLSPLTPLTVVTLLQASNPAAVVVGRLLQAATNYHNGHTGQLSAITVFLFGGSLARIFTSIQETGD PLMAGTFVVSSLCNGLIATQLLFYWNAKPPHKQKKAQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_004870
ORF Size:	741 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_004870.2</u>
RefSeq Size:	1634 bp
RefSeq ORF:	744 bp
Locus ID:	9526
UniProt ID:	<u>O75352</u>
Cytogenetics:	17p13.1
MW:	26.7 kDa