

## Product datasheet for **RC200576L2V**

### **HNRPUL1 (HNRNPUL1) (NM\_007040) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	HNRPUL1 (HNRNPUL1) (NM_007040) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HNRPUL1
Synonyms:	E1B-AP5; E1BAP5; HNRPUL1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_007040
ORF Size:	2568 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200576).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007040.2</a>
RefSeq Size:	3892 bp
RefSeq ORF:	2571 bp
Locus ID:	11100
UniProt ID:	<a href="#">Q9BUJ2</a>
Cytogenetics:	19q13.2
Domains:	SAP, SPRY
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 95.8 kDa

**Gene Summary:** This gene encodes a nuclear RNA-binding protein of the heterogeneous nuclear ribonucleoprotein (hnRNP) family. This protein binds specifically to adenovirus early-1B-55kDa oncoprotein. It may play an important role in nucleocytoplasmic RNA transport, and its function is modulated by early-1B-55kDa in adenovirus-infected cells. [provided by RefSeq, Mar 2016]