

Product datasheet for **RC200586L1V**

FSTL1 (NM_007085) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FSTL1 (NM_007085) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FSTL1
Synonyms:	FRP; FSL1; MIR198; OCC-1; OCC1; tsc36
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_007085
ORF Size:	924 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200586).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007085.3
RefSeq Size:	3840 bp
RefSeq ORF:	927 bp
Locus ID:	11167
UniProt ID:	Q12841
Cytogenetics:	3q13.33
Domains:	EFh, kazal, FOLN
Protein Families:	Secreted Protein



[View online »](#)

MW: 35 kDa

Gene Summary: This gene encodes a protein with similarity to follistatin, an activin-binding protein. It contains an FS module, a follistatin-like sequence containing 10 conserved cysteine residues. This gene product is thought to be an autoantigen associated with rheumatoid arthritis. [provided by RefSeq, Jul 2008]