

Product datasheet for **RC200587L1V**

RAMP1 (NM_005855) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | RAMP1 (NM_005855) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | RAMP1 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_005855 |
| ORF Size: | 444 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC200587). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_005855.2 |
| RefSeq Size: | 922 bp |
| RefSeq ORF: | 447 bp |
| Locus ID: | 10267 |
| UniProt ID: | O60894 |
| Cytogenetics: | 2q37.3 |
| Domains: | RAMP |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | Vascular smooth muscle contraction |



[View online »](#)

MW: 17 kDa

Gene Summary: The protein encoded by this gene is a member of the RAMP family of single-transmembrane-domain proteins, called receptor (calcitonin) activity modifying proteins (RAMPs). RAMPs are type I transmembrane proteins with an extracellular N terminus and a cytoplasmic C terminus. RAMPs are required to transport calcitonin-receptor-like receptor (CRLR) to the plasma membrane. CRLR, a receptor with seven transmembrane domains, can function as either a calcitonin-gene-related peptide (CGRP) receptor or an adrenomedullin receptor, depending on which members of the RAMP family are expressed. In the presence of this (RAMP1) protein, CRLR functions as a CGRP receptor. The RAMP1 protein is involved in the terminal glycosylation, maturation, and presentation of the CGRP receptor to the cell surface. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015]