

Product datasheet for **RC200600L1V**

VGLL1 (NM_016267) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	VGLL1 (NM_016267) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VGLL1
Synonyms:	TDU; VGL1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016267
ORF Size:	774 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200600).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016267.2
RefSeq Size:	1242 bp
RefSeq ORF:	777 bp
Locus ID:	51442
UniProt ID:	Q99990
Cytogenetics:	Xq26.3
Domains:	TDU
Protein Families:	Transcription Factors



[View online »](#)

MW: 28.7 kDa

Gene Summary: The protein encoded by this gene binds proteins of the TEA domain family of transcription factors (TEFs) through the Vg (vestigial) homology region found in its N-terminus. It may thus function as a specific coactivator for the mammalian TEFs. [provided by RefSeq, Sep 2009]