

Product datasheet for **RC200602A1V**

Human TRAPPC3 (NM_014408) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TRAPPC3 (NM_014408) AAV Particle
Tag: Myc-DDK
Symbol: TRAPPC3
Synonyms: BET3
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200602 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCGAGGCAGGCGAACCGTGGCACCAGAGCAAGAAAATGAGCTCTGAGCTCTTCACCTGACCTATG
 GTGCCCTGGTCACCCAGCTATGTAAGGACTATGAAAATGATGAAGATGTGAATAAACAGCTGGACAAAAT
 GGGCTTTAACATTGGAGTCCGGCTGATTGAAGATTTCTGGCTCGGTCAAATGTTGGGAGGTGCCATGAC
 TTTCGGGAAACTGCGGATGTCAATGCCAAGGTGGCGTTCAAGATGTACTTGGGCATCACTCCAAGCATT
 CTAATTGGAGCCAGCTGGTGATGAATTCCTCCCTCATTTTGAAAATAACCCCTTGGTGGACTTTGTGGA
 ACTTCCTGATAACCACTCATCCCTTATTTATTCCAATCTCTTGTGTGGGTGTTGCGGGGAGCTTTGGAG
 ATGGTCCAGATGGCTGTGGAGGCCAAGTTGTCCAGGACACCCTGAAAGGAGACGGTGTGACAGAAATCC
 GGATGAGATTCATCAGGCGGATTGAGGACAATCTCCAGCTGGAGAGGAA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200602 protein sequence
 Red=Cloning site Green=Tags(s)

MSRQANRGTESKKMSELFLLTYGALVTQLCKDYENDEVDNKLQDKMGFNIGVRLIEDFLARSNVGRCHD
 FRETADVIAKVAFKMYLGITPSITNWPAGDEFSLILENNPLVDFVELPDNHSSLIYSNLLCGVLRGALE
 MVQMAVEAKFVQDTLKGDTVTEIRMRFIRRIEDNLPAGEE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_014408
ORF Size:	540 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014408.3</u>
RefSeq Size:	1333 bp
RefSeq ORF:	543 bp
Locus ID:	27095
UniProt ID:	<u>O43617</u>
Cytogenetics:	1p34.3
MW:	20.3 kDa